

## Taq DNA polymerase Unis

(Catalogue number T 037, T 038, T 039)

rev. 04/2025

### Description

Taq DNA polymerase Unis is an alternative to Taq DNA polymerase (Cat. No. T032-T034). Difference between these two products is in difference in enzyme storage buffer which enhances stability and efficiency of the enzyme. Occasional amplification of nonspecific fragments can be prevented by enhanced dilution of the enzyme.

### Taq DNA polymerase

*Taq DNA polymerase is a thermostable enzyme isolated from Thermus aquaticus. The enzyme catalyzes synthesis of complementary DNA strand in the 5'->3' direction and also possesses a 5'->3' exonuclease activity. During amplification of DNA fragments, Taq polymerase adds at 3' end an adenosine overhang. This can be utilized for cloning of PCR-generated DNA fragments. Advantage of the enzyme is its high processivity (amplification of 1000 base pairs takes < 1 min). Disadvantage of the enzyme is that it lacks a 3'->5' exonuclease proofreading activity and this accounts for high error rate [about 1 error to 10<sup>5</sup> - 10<sup>6</sup> base pairs (bps)]. The major usage of the enzyme is in diagnostic analysis for amplification of DNA fragments up to 5000 bps.*

### Technical data

#### Components and packaging

- Taq DNA polymerase Unis is supplied at a concentration 5 U/μl. Basic packaging is 1 tube with 500 U/100 μl (T037), 5 tubes with 500 U/100 μl (T038) or 10 tubes with 500 U/100 μl (T039).
- Each tube of Taq DNA polymerase Unis is accompanied by a tube with 10x concentrated react buffer containing MgCl<sub>2</sub> (1.5 ml). If different concentration of MgCl<sub>2</sub> is required, a tube with 10x concentrated reaction buffer without MgCl<sub>2</sub> (1.5 ml) and a tube with 25 mM MgCl<sub>2</sub> (0.5 ml) should be ordered (Cat. No. T035).

#### Storage

- At temperature -20°C ± 5°C. Material can be repeatedly defrosted.

#### Composition

- Storage buffer for Taq DNA polymerase Unis: 20 mM HEPES (pH 7.9 at 25°C), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5 mM PMSF, stabilizers, 50% glycerol.
- 10x reaction buffer: 100 mM Tris-HCl, pH 8.8 (at 25°C), 500 mM KCl, 1% Triton X-100, 15 mM MgCl<sub>2</sub>.

#### Activity

- One unit is defined as the amount of enzyme, which catalyzes incorporation of 10 nmol dNTPs within 30 min at 72°C into trichloroacetic acid precipitable material. Reaction conditions are as follow: 10 mM Tris-HCl (pH 8.8 at 25°C), 50 mM KCl, 0.1% Triton X-100, 1.5 mM MgCl<sub>2</sub>, 200 μM dATP, dCTP, dGTP and [α-<sup>32</sup>P]dTTP, 50 μg/ml denatured cDNA, 0.5 μM primer and 0.2 – 0.5 U of enzyme in the volume of 50 μl.

#### Purity and quality control

- Purity of Taq DNA polymerase Unis is tested by SDS-PAGE. After staining with Coomassie blue, enzyme migrates as the major band of 94 kDa. Material is nuclease free.
- Each batch of Taq DNA polymerase Unis is tested for its ability to amplify DNA fragment from mammalian genomic DNA by PCR. The results are verified by electrophoresis in agarose gel in the presence of ethidium bromide; only DNA band of the expected size is present

Cat. No.	Product name and specification	Amount
T037	Taq DNA polymerase Unis	500U
T038	Taq DNA polymerase Unis	5x 500U
T039	Taq DNA polymerase Unis	10x 500U
T035	10x conc. Taq buffer without MgCl <sub>2</sub> +MgCl <sub>2</sub>	1.5 ml + 0.5 ml



# Taq DNA polymerase Unis

(Catalogue number T 037, T 038, T 039)

## Protocol

### Basic protocol

The protocol described below can be used for routine PCR. However, in some cases reaction conditions must be optimized, mainly annealing temperature and concentration of MgCl<sub>2</sub>.

1. In thin-wall test tubes the following components are mixed<sup>1</sup>:

	PCR in 50 µl	Final concentration
10x reaction buffer with MgCl <sub>2</sub> <sup>2</sup>	5 µl	1x react. buffer with 1.5 mM MgCl <sub>2</sub>
PCR dNTP mix (10 mM each) (Cat. No. P041)	1 µl	0.2 mM dNTP each
5' primer (50 µM)	0.5 µl	0.5 µM
3' primer (50 µM)	0.5 µl	0.5 µM
Taq DNA polymerase Unis (5U/µl)	0.5 µl	2.5 U (0.05 U/µl)
Template DNA (1 ng/µl - 1 µg/µl)	1 µl	0.02 ng/µl – 0.02 µg/µl
PCR H <sub>2</sub> O (Cat. No. P042)	41.5 µl	

<sup>1</sup> When more DNA samples are tested with the same primers, it is convenient to prepare so called Master Mix in which the volumes of individual components are multiplied by the number of DNA samples tested. Then, 49 µl aliquots of the Master Mix are distributed into the PCR tubes, followed by addition of tested DNA (1 µl) into each tube.

<sup>2</sup> If unsatisfactory results are obtained, PCRs with different concentrations of MgCl<sub>2</sub> should be performed (see below).

2. Samples are homogenized and spun down.

3. If cycler without heating lid is used, PCR oil 25 µl (Cat. No. P043) is added.

4. Thermal cycler is programmed according to the manufacturer's instructions. A typical cycling program is:

	Temperature	Time	Number of cycles
Initial denaturation	94°C	3 min	1
Denaturation	94°C	30 s	25-35
Annealing of primers	55-68°C <sup>1</sup>	30 s	
Extension	72°C	1 min per 1 kb	
Final extension	72°C	10 min	1
Cooling	4°C		

<sup>1</sup> Should be determined experimentally; usually 5°C below melting temperature of the primers.

5. After PCR, samples are mixed with loading buffer (Cat. No. P048, P062, P066 or P065) and analyzed by electrophoresis in agarose gel in the presence of ethidium bromide. Alternatively, the samples can be stored at -20°C.

### Optimization of MgCl<sub>2</sub> concentration

MgCl<sub>2</sub> at a 1.5 mM final concentration is suitable for most PCRs. However, if amplification of nonspecific DNA fragments is observed, optimal Mg<sup>2+</sup> concentration for given PCR should be determined. To this end 10x reaction buffer without MgCl<sub>2</sub> and 25 mM MgCl<sub>2</sub> can be ordered (Cat. No. T035).

1. Preparation of Master Mix **without MgCl<sub>2</sub>** by mixing the following components:

10x reaction buffer without MgCl <sub>2</sub>	40 µl
PCR dNTP mix (10 mM each)	8 µl
5' primer (50 µM)	4 µl
3' primer (50 µM)	4 µl
Taq DNA polymerase Unis (5U/µl)	4 µl
Template DNA (1 ng/µl - 1 µg/µl)	8 µl
PCR H <sub>2</sub> O	268 µl
<b>Total volume</b>	<b>336 µl</b>

2. Master mix is thoroughly mixed, centrifuged briefly and 42  $\mu\text{l}$  aliquots are distributed into 7 PCR tubes.
3.  $\text{MgCl}_2$  and PCR  $\text{H}_2\text{O}$  is added into PCR Master mixes as follows:

<b>Tube No.</b>	<b>25 mM <math>\text{MgCl}_2</math></b>	<b>PCR <math>\text{H}_2\text{O}</math></b>	<b>Final <math>\text{MgCl}_2</math> concentration</b>
1	1 $\mu\text{l}$	7 $\mu\text{l}$	0.5 mM
2	2 $\mu\text{l}$	6 $\mu\text{l}$	1.0 mM
3	3 $\mu\text{l}$	5 $\mu\text{l}$	1.5 mM
4	4 $\mu\text{l}$	4 $\mu\text{l}$	2.0 mM
5	5 $\mu\text{l}$	3 $\mu\text{l}$	2.5 mM
6	6 $\mu\text{l}$	2 $\mu\text{l}$	3.0 mM
7	8 $\mu\text{l}$	0 $\mu\text{l}$	4.0 mM

4. PCR is performed as above and the samples are analyzed by electrophoresis in agarose gel in the presence of ethidium bromide. Optimal concentration of  $\text{MgCl}_2$  for given PCR is thus determined.