

PCR ethidium bromide

IN QUALITY SUITABLE FOR THE ANALYSIS OF PCR PRODUCTS
(Cat. No. P046, P047)

rev. 01/2022

Description

It is a filtered solution of ethidium bromide (10 mg/ml) in PCR H₂O.

Technical data

Storage

- Store at temperature $4 \pm 3^{\circ}\text{C}$.

Packaging

- 1 ml in plastic tube with screw cap.

Quality control

- Each batch of PCR ethidium bromide is tested for its ability to visualize DNA fragments by electrophoresis in PCR agarose gel (Cat. No. P045). The final concentration of ethidium bromide in gel is usually 0.5 $\mu\text{g/ml}$.

Cat. No.	Product name and specification	Amount
P046	PCR ethidium bromide	1ml
P047	PCR ethidium bromide	5x 1ml

