

**CERTIFICATE OF ANALYSIS**

**Product:** qPCR 2x SYBR Master Mix

**Catalog No:** P551, P552, P553, P553xl

**Lot No:** P551042024

**Date of Expiry:** 04/2024

**Composition:** qPCR 2x SYBR Master Mix contains: 20 mM Tris-HCl, pH 8.8 (at 25°C), 100 mM KCl, 0.2% Triton X-100, 3 mM MgCl<sub>2</sub>, 400 μM dATP, 400 μM dCTP, 400 μM dGTP, 400 μM dTTP, 50 U/ml Taq DNA polymerase, monoclonal antibody anti-Taq (38 nM), SYBR Green I, stabilizers and additives.

**Supplied with:** PCR Ultra H<sub>2</sub>O (Cat. No. P040).

**Storage temperature:** For short terms (days) at 4°C ± 3°C.  
For long terms at -20 ± 5°C. Material can be repeatedly defrosted.

**Functional Test:** The lot has been tested for the ability to amplify a fragment of genomic DNA using the following conditions:

**Test conditions:**

Volume*	Reagent	Final concentration
12.5 μl	qPCR 2x SYBR Master Mix	1x qPCR 2x SYBR Master Mix
0.5 μl	Forward primer	50 μM 5' primer, 5'-TTGGAGAGTTTGATCCTGGCTC-3'
0.5 μl	Reverse primer	50 μM 3' primer 5'-TTGGAGAGTTTGATCCTGGCTC-3'
1 μl	Template DNA	containing 10 ng of E-coli DNA
10.5 μl	PCR Ultra H <sub>2</sub> O	(to a final volume 25 μl)

**Cycling conditions on LightCycler 96, ROCHE:**

	Temperature	Time	Number of cycles
Initial denaturation	94°C	3 min	1
Denaturation	94°C	10 s	35
Annealing of primers	55°C	10 s	
Extension	72°C	10 s	
Final extension	72°C	7 min	1
Cooling	22°C		

**Result:** As expected, Ct values were between 18 - 20, whereas in the absence of DNA Ct >35. Furthermore, electrophoresis of the PCR product on agarose gel revealed one band of 173 bp.

**FOR RESEARCH USE****APPROVED DATE:** 28.10.2021

Manager: Hana Těšitelová